

CHIRANA

CONNECT



Reliability + Precision + Comfort = Chirana.

CHIRANA CONNECT

Designed for maximum efficiency, patient comfort, and seamless workflow integration, this dental unit combines modern aesthetics with reliable performance.



Slovak manufacturer

We are a proud Slovak manufacturer with decades of expertise in dental technology. Our newest dental unit, the CHIRANA CONNECT, is designed and manufactured in Slovakia to the highest European standards.

Our commitment doesn't end with delivery

We provide exceptional after-sales support to keep your practice running smoothly. With fast response times and expert service, we ensure your CHIRANA CONNECT performs flawlessly for years."



Designed for maximum efficiency, patient comfort, and seamless workflow integration.



Dentist Table

Intuitive work space for you

Compact Size. Big Performance

Our new compact dentist table delivers full functionality in a space-saving design

Key Features

- **Five Instrument Holders** for enhanced treatment versatility
- **Continental** (Upper) or **Traditional** (Lower) Hose Delivery options
- **Preset Settings** for up to 4 users
- **Water Regulation** – single regulator or independent regulators for each instrument
- **Electronic Water** and **Air Regulation** for precision control
- **Manual** or **Pneumatic Brake System** for secure positioning
- **Glass Touchscreen Keyboard** – intuitive control, easy cleaning, and superior hygiene



Assistant table

Closer to the Patient



Key Features

- **Three Instrument Holders** for versatile treatment options
- **Touch Glass Keyboard** for intuitive control
- **Simplified Arm** for maximum space efficiency
- **Liquid Detection Sensors** for hygiene and distilled water monitoring
- **Connected Programming directly** from the dentist table
- **Integrated Functions** – program cup filling and bowl flush
- **Pantographic Arm** for optimal reach and flexibility



Dental chair

Designed for Pure Comfort



Key Features

- **Trendelenburg Position** for optimal patient comfort and clinical access
- **2D or 3D Headrest** for precise head and neck support
- **Ergonomic Left Armrest** (standard) and **Foldable Right Armrest** for easy access
- **Synchronized Movement** for smooth, coordinated adjustments
- **Programmable Positions** for up to 4 users
- **Triple Control Options** – dentist table, assistant table, or multifunctional foot controller
- **Upholstery Choices** – standard, stitched, or premium for maximum comfort
- **Wide Color Selection** to match any practice style



Spitton block

Accessible Maintenance
Anytime



Key Features

- **Ceramic Bowl** – easy to remove and clean
- **Programmable Functions** – cup filling and bowl flush
- **Removable Cover** – easy access to service space
- **Compact Design** – saves space without compromising functionality
- **Versatile Suction Options** – compatible with wet or dry ejectors (DURR, CATANI, Chirana, etc.)
- **Dual-Side Access** – one side dedicated to routine day-to-day use, the other optimized for full-service maintenance



Connect options

Closer to the Patient




Key Features

- **External Media Connections** – compatible with any tabletop solution
- **Multifunctional Foot Controller** – wired or wireless operation
- **Integrated Programming** – control across the entire dental unit
- **Versatile Platform** – allows customization to meet various requirements
- **Complete Your Dental Unit** – fill it with our range of Dental instruments for a fully equipped solution



No tools required

"Open your dental unit without any tools. Simplify maintenance like never before."



CHIRANA
Medical STARÁ
TURÁ



**connect
yourselft
with
comfort**

"Discover more details on
our website."

www.chirana.eu